


Fetal Echocardiography

Pediatric Echocardiography Reporting Module

ACV 8.0


doc v1

The pediatric echocardiography module contains much of the same content as the adult echocardiography module but is enriched to support fetal, infant, and child studies.






When you begin reporting in ASCEND, selections made in the *Startup* screen build a tab set uniquely tailored to that study. ASCEND is designed to support a single fetal study or multiple pregnancies. For this example we will choose two fetuses.




The *Fetal* tab lists the cardiac structures within their finding group. Every finding hosts a *Free text* note icon as a quick and easy way to insert a dictated comment into the report.





The *Findings per fetus* tab gives you a way to identify and describe each fetus individually. By selecting the *Fetal description* data entry form, detailed measurements, heart rhythms and overall impressions can be reported for each fetus.





www.ascendhit.com

HEALTH INFORMATION TECHNOLOGY

Phone (Toll Free): 844-413-2610

Email: information@ascendhit.com